

MELTING OF EVA UNDER CO₂ PRESSURE FOR SOLAR PANEL RECYCLING

INTRODUCTION


Apart from the active photovoltaic (PV) material, solar cells are made of several components with separate functions. The PV material is typically encapsulated in a transparent material, usually EVA polymer, to protect it from deterioration.

Recycling such layered material is challenging. A process to separate PV material from EVA has been proposed by researchers. It consists in inducing delamination by applying high pressure of CO₂. In this study, the effect of CO₂ pressure on melting properties of EVA was studied using a Calvet Pro instrument.

EXPERIMENT

- Sample: EVA-28 (10mg)
- Cell: High pressure cell connected to a high pressure syringe pump
- Instrument: Calvet Pro (formerly Sensys Evo)
- Condition: Heating to 10°C/min
- Atmosphere: CO₂ (1bar, 30bar, 60bar, 100bar, 150bar)

RESULTS



Upon heating, a melting endotherm is measured above 65°C. The onset temperature and the width of the melting peak are modified by the CO₂ pressure.


Pressure induces a decrease in the melting range. This is explained by its impact on the polymer crystallinity and the presence of CO₂.

This test is useful to evaluate the effect of CO₂ pressure on the structure and transformation of EVA and to select the appropriate process conditions for solar panel recycling.

Axel Briand et al., CO₂ absorption into a polymer within a multilayer structure: The case of poly(ethylene-co-vinyl acetate) in photovoltaic modules, The Journal of Supercritical Fluids, Volume 179, 2022, 105380, ISSN 0896-8446

INSTRUMENT

CALVET PRO



- HIGHEST HEAT MEASUREMENT ACCURACY**
Calvet 3D sensor based on thermocouples with Joule effect calibration
- HIGHEST MASS VARIATION ACCURACY**
with its optional Hang-down Symmetrical Beam Balance
- SUB-AMBIENT TO HIGH TEMPERATURE OPERATIONS**
with solutions from -120 to 830 °C
- ISOTHERMAL OR TEMPERATURE SCANNING MODES**
for increased flexibility and replication of real life conditions
- CONVENIENT INTERCHANGEABLE CRUCIBLES AND CELLS**
to perform even the most demanding experiments using one instrument :
 - high pressure (up to 500 bar) and high vacuum
 - pressure resistance, measurement or control
 - packed bed reactor experiments
- EXTERNAL COUPLING CAPABILITY**

Switzerland – France – China – United States – India – Hong Kong - www.setaram.com - sales@setaram.com

SETARAM is a registered trademark of KEP Technologies Group

