FLEXIWET



HIGH VERSATILITY AND ACCURACY WET GAS GENERATOR

for multiple laboratory or industrial applications

- Accurate wet conditions control
- From simple setpoint to sophisticated programming
- Operation with various types of gases

PLUG&PLAY, SIMPLICITY, EASY SETTINGS

- Easily transportable
- Simple connection system to any laboratory instruments or climate chambers
- Simple and convenient user interface

PERFORMANCE (GAS)		
Pre-calibrated for various gases		Air, helium, nitrogen, carbon dioxide, argon
Flow Rate	FLEXI WET 50	3 to 50 ml/min
	FLEXI WET 200	10 to 200 ml/min
Heated Transfer line		Ambient to 100°C
PERFORMANCE (HUMIDITY)		
Gas humidity	Ambient to 50 °C	0% RH°; 5-95% RH
	50 to 70 °C	0% RH³; 5-90% RH
Humidity profile generation		Constant RH, steps, ramps
Autonomy		> 1000 hours ^b
Accuracy		+/- 0.8% RH
Accuracy		+/- 0.1 °C
Stability		+/- 0.3% RH
External RH probe		Optional
GENERAL		
Size (W x D x H)		420 x 530 x 350 mm 16.5 x 20.9 x 13.8 inch
Weight		22 kg 49 lb
Power supply		110 / 230 V 50/60 Hz

^aoperations with dry gas: isolation of the saturator using an automatic valve switch, ^bat 70°C, 90% RH and 20 ml/min